

A16450

Leader in Biomolecular Solutions for Life Science



EPX Rabbit pAb

Catalog No.: A16450

Basic Information

Observed MW

81kDa

Calculated MW

81kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene is a member of the peroxidase gene family and is expressed in eosinophils. The encoded preproprotein is proteolytically processed into covalently attached heavy and light chains to form the mature enzyme, which functions as an oxidant. The enzyme is released at sites of parasitic infection or allergen stimulation to mediate lysis of protozoa or parasitic worms. The gene is found in a gene cluster with other peroxidase genes on chromosome 17. Mutations in this gene result in eosinophil peroxidase deficiency.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8288

Swiss Prot

P11678

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 400-500 of human EPX (NP_000493.1).

Synonyms

EPO; EPP; EPXD; EPX-PEN; EPX

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using EPX Rabbit pAb (A16450) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.