

A16515

Leader in Biomolecular Solutions for Life Science



BRPF3 Rabbit pAb

Catalog No.: A16515

Basic Information

Observed MW

Refer to figures

Calculated MW

136kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Rat

Background

Predicted to enable histone binding activity. Involved in histone H3-K14 acetylation and positive regulation of DNA replication. Part of MOZ/MORF histone acetyltransferase complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

27154

Swiss Prot

Q9ULD4

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human BRPF3 (NP_056510.2).

Synonyms

BRPF3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.