

A16580

Leader in Biomolecular Solutions for Life Science



FAM234A Rabbit pAb

Catalog No.: A16580

Basic Information

Observed MW

67kDa

Calculated MW

60kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Located in cell surface.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

83986

Swiss Prot

Q9H0X4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 70-400 of human FAM234A (NP_114428.1).

Synonyms

gs19; ITFG3; C16orf9; FAM234A

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using FAM234A Rabbit pAb (A16580) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 3s.