

CNN1 Rabbit pAb

Catalog No.	A16638	Category	Polyclonal Antibodies
Applications	WB	Observed MW	33kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	31kDa/33kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-180 of human CNN1 (NP_001290.2).
Gene ID	1264
Swiss prot	P51911
Synonyms	CNN1;HEL-S-14;SMCC;Sm-Calp


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 - 1:2000



Western blot analysis of extracts of various cell lines, using CNN1 antibody (A16638) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.