

A16754

Leader in Biomolecular Solutions for Life Science



## XIAP Rabbit pAb

Catalog No.: A16754

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

57kDa

#### Category

Polyclonal Antibody

#### Applications

WB, IF/ICC

#### Cross-Reactivity

Human, Mouse, Rat

### Background

This gene encodes a protein that belongs to a family of apoptotic suppressor proteins. Members of this family share a conserved motif termed, baculovirus IAP repeat, which is necessary for their anti-apoptotic function. This protein functions through binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2 and inhibits apoptosis induced by menadione, a potent inducer of free radicals, and interleukin 1-beta converting enzyme. This protein also inhibits at least two members of the caspase family of cell-death proteases, caspase-3 and caspase-7. Mutations in this gene are the cause of X-linked lymphoproliferative syndrome. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 2 and 11.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:100

### Immunogen Information

#### Gene ID

331

#### Swiss Prot

P98170

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human XIAP (NP\_001158.2).

#### Synonyms

API3; ILP1; MIHA; XLP2; BIRC4; IAP-3; hIAP3; hIAP-3; XIAP

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.