

NPAT Rabbit pAb

Catalog No.: A16931

Basic Information

Observed MW

Refer to figures

Calculated MW

154kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables protein C-terminus binding activity; transcription coactivator activity; and transcription corepressor activity. Involved in positive regulation of transcription by RNA polymerase II and regulation of transcription involved in G1/S transition of mitotic cell cycle. Located in Cajal body; Gemini of coiled bodies; and cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4863

Swiss Prot

Q14207

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1128-1427 of human NPAT (NP_002510.2).

Synonyms

E14; p220; E14/NPAT; NPAT

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.