

A17209

Leader in Biomolecular Solutions for Life Science



Bcl9 Rabbit pAb

Catalog No.: A17209

Basic Information

Observed MW

Refer to figures

Calculated MW

149kDa

Category

Polyclonal Antibody

Applications

WB, IHC-P, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables beta-catenin binding activity. Acts upstream of or within several processes, including myotube differentiation involved in skeletal muscle regeneration; skeletal muscle cell differentiation; and somatic stem cell population maintenance. Located in cis-Golgi network and nucleus. Is expressed in several structures, including brain; eye; head surface ectoderm; metanephros; and thymus primordium. Orthologous to human BCL9 (BCL9 transcription coactivator).

Recommended Dilutions

WB	1:500 - 1:2000
IHC-P	1:50 - 1:200
IF/ICC	1:50 - 1:200

Immunogen Information

Gene ID

77578

Swiss Prot

Q9D219

Immunogen

Recombinant protein of mouse Bcl9.

Synonyms

Gm130; 2610202E01Rik; 8030475K17Rik; A330041G23Rik; Bcl9

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Immunofluorescence analysis of U2OS cells using Bcl9 Rabbit pAb (A17209) at dilution of 100 (40x lens). Blue: DAPI for nuclear staining.