

A17348

Leader in Biomolecular Solutions for Life Science



ZDHC8 Rabbit pAb

Catalog No.: A17348

Basic Information

Observed MW

81kDa

Calculated MW

81kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

This gene encodes a four transmembrane protein that is a member of the zinc finger DHHC domain-containing protein family. The encoded protein may function as a palmitoyltransferase. Defects in this gene may be associated with a susceptibility to schizophrenia. Alternate splicing of this gene results in multiple transcript variants. A pseudogene of this gene is found on chromosome 22.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

29801

Swiss Prot

Q9ULC8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 240-380 of human ZDHC8 (NP_037505.1).

Synonyms

DHHC8; ZNF378; ZDHC1; ZDHC8

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data

