

A17372

Leader in Biomolecular Solutions for Life Science



## HERC4 Rabbit pAb

Catalog No.: A17372

1 Publications

### Basic Information

#### Observed MW

119kDa

#### Calculated MW

119kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse, Rat

### Background

HERC4 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

26091

#### Swiss Prot

Q5GLZ8

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human HERC4 (NP\_056416.2).

#### Synonyms

HERC4

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.