

AZI2 Rabbit pAb

Catalog No.: A17756

Basic Information

Observed MW

Refer to figures

Calculated MW

45kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

AZI2, or NAP1, contributes to the activation of NF κ B (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK (TBK1; MIM 604834) (Fujita et al., 2003 [PubMed 14560022]).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

64343

Swiss Prot

Q9H6S1

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 100-200 of human AZI2 (NP_071906.1).

Synonyms

AZ2; NAP1; TILP; AZI2

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.