

A17861

Leader in Biomolecular Solutions for Life Science



FILIP1L Rabbit pAb

Catalog No.: A17861

Basic Information

Observed MW

110kDa

Calculated MW

130kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to be located in cytoplasm; membrane; and nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

11259

Swiss Prot

Q4L180

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FILIP1L (NP_878913.2).

Synonyms

DOC1; DOC-1; GIP90; GIP130; FILIP1L

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using FILIP1L Rabbit pAb (A17861) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.