

A18153

Leader in Biomolecular Solutions for Life Science



Ajuba Rabbit pAb

Catalog No.: A18153

Basic Information

Observed MW

55kDa

Calculated MW

57kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables alpha-catenin binding activity and transcription corepressor activity. Involved in several processes, including negative regulation of hippo signaling; positive regulation of gene silencing by miRNA; and regulation of cellular response to hypoxia. Acts upstream of or within gene silencing by miRNA and positive regulation of protein-containing complex assembly. Located in several cellular components, including Golgi apparatus; P-body; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

84962

Swiss Prot

Q96IF1

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Ajuba (NP_116265.1).

Synonyms

JUB; Ajuba

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.