

A18271

Leader in Biomolecular Solutions for Life Science



## OSR2 Rabbit pAb

Catalog No.: A18271

### Basic Information

**Observed MW**

38kDa

**Calculated MW**

36kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human

### Background

OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).

### Recommended Dilutions

**WB** 1:500 - 1:2000

### Immunogen Information

**Gene ID**

116039

**Swiss Prot**

Q8N2R0

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human OSR2 (NP\_001135934.1).

**Synonyms**

OSR2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of lysates from 293T cells, using OSR2 Rabbit pAb (A18271) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 180s.