

ANKRD26 Rabbit pAb

Catalog No.: A18422

Basic Information

Observed MW

Refer to figures

Calculated MW

196kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse

Background

This gene encodes a protein containing N-terminal ankyrin repeats which function in protein-protein interactions. Mutations in this gene are associated with autosomal dominant thrombocytopenia-2. Pseudogenes of this gene are found on chromosome 7, 10, 13 and 16. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

22852

Swiss Prot

Q9UP58

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 320-440 of human ANKRD26 (NP_055730.2).

Synonyms

THC2; bA145E8.1; ANKRD26

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.