

A18442

Leader in Biomolecular Solutions for Life Science



## BRI3 Rabbit pAb

Catalog No.: A18442

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

14kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

Enables identical protein binding activity. Predicted to be located in azurophil granule membrane and plasma membrane.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

25798

#### Swiss Prot

O95415

#### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BRI3 (NP\_056194.1).

#### Synonyms

I3; BRI3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.