

A18475

Leader in Biomolecular Solutions for Life Science



ZNF302 Rabbit pAb

Catalog No.: A18475

Basic Information

Observed MW

55kDa

Calculated MW

55kDa

Category

Antibody

Applications

WB,ELISA

Cross-Reactivity

Rat

Background

This gene encodes a member of the zinc-finger protein family. The encoded protein contains seven C2H2-type zinc fingers and a KRAB domain, but its function has yet to be determined. Alternatively spliced transcript variants have been described.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

55900

Swiss Prot

Q9NR11

Immunogen

Recombinant Protein corresponding to a sequence within amino acids 90-200 of human ZNF302 (NP_060913.2).

Synonyms

HSD16; MST154; ZNF327; MSTP154; ZNF135L; ZNF140L; ZNF302

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from Rat testis using ZNF302 Rabbit pAb(A18475) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.