

A18516

Leader in Biomolecular Solutions for Life Science



## FBF1 Rabbit pAb

Catalog No.: A18516

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

125kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Rat

### Background

Involved in apical junction assembly and establishment of epithelial cell polarity. Acts upstream of or within cilium assembly. Located in centriole; centrosome; and spindle pole. Part of ciliary transition fiber. Colocalizes with apical junction complex and keratin filament.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

85302

#### Swiss Prot

Q8TES7

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human FBF1 (Q8TES7).

#### Synonyms

Alb; FBF-1; FBF1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.