

A18546

Leader in Biomolecular Solutions for Life Science



## PLEKHA7 Rabbit pAb

Catalog No.: A18546

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

127kDa

#### Category

Polyclonal Antibody

#### Applications

WB,IF/ICC,ELISA

#### Cross-Reactivity

Human,Mouse,Rat

### Background

Enables delta-catenin binding activity. Involved in epithelial cell-cell adhesion; pore complex assembly; and zonula adherens maintenance. Located in several cellular components, including centrosome; nucleoplasm; and zonula adherens. Part of pore complex.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

### Immunogen Information

#### Gene ID

144100

#### Swiss Prot

Q6IQ23

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PLEKHA7 (NP\_778228.3).

#### Synonyms

PLEKHA7

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.