

A18656

Leader in Biomolecular Solutions for Life Science



BASP1 Rabbit pAb

Catalog No.: A18656

Basic Information

Observed MW

Refer to figures

Calculated MW

23kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a membrane bound protein with several transient phosphorylation sites and PEST motifs. Conservation of proteins with PEST sequences among different species supports their functional significance. PEST sequences typically occur in proteins with high turnover rates. Immunological characteristics of this protein are species specific. This protein also undergoes N-terminal myristoylation. Alternative splicing results in multiple transcript variants that encode the same protein.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

10409

Swiss Prot

P80723

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human BASP1 (NP_006308.3).

Synonyms

CAP23; NAP22; CAP-23; NAP-22; BASP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.