

A18822

Leader in Biomolecular Solutions for Life Science



Methyl-SUPT5H/SPT5-H418 Rabbit pAb

Catalog No.: A18822

Basic Information

Observed MW

Refer to figures

Calculated MW

121kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables enzyme binding activity and protein heterodimerization activity. Involved in positive regulation of macroautophagy; regulation of RNA metabolic process; and transcription elongation from RNA polymerase II promoter. Located in nucleoplasm. Part of DSIF complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

6829

Swiss Prot

O00267

Immunogen

A synthetic methylated peptide around H418 of human SUPT5H/SPT5 (NP_003160.2).

Synonyms

SPT5; SPT5H; Tat-CT1; Methyl-SUPT5H/SPT5-H418

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.