

A18864

Leader in Biomolecular Solutions for Life Science



TMEM187 Rabbit pAb

Catalog No.: A18864

Basic Information

Observed MW

Refer to figures

Calculated MW

29kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene consists of two exons and encodes a multi-pass membrane protein. An alternatively spliced transcript variant encoding the same protein has been found, but its biological validity is not determined.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

8269

Swiss Prot

Q14656

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMEM187 (NP_003483.1).

Synonyms

ITBA1; CXorf12; DXS9878E; TMEM187

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.