

A18867

Leader in Biomolecular Solutions for Life Science



HIVp24 Rabbit pAb

Catalog No.: A18867

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Species independent

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-231 of human immunodeficiency virus 1 immunodeficiency virus 1 HIVp24 (NP_579880.1).

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.