

A1887

Leader in Biomolecular Solutions for Life Science



DNA Ligase III/LIG3 Rabbit pAb

Catalog No.: A1887

3 Publications

Basic Information

Observed MW

113kDa

Calculated MW

113kDa

Category

Polyclonal Antibody

Applications

WB, IP, ELISA

Cross-Reactivity

Human

Background

This gene is a member of the DNA ligase family. Each member of this family encodes a protein that catalyzes the joining of DNA ends but they each have a distinct role in DNA metabolism. The protein encoded by this gene is involved in excision repair and is located in both the mitochondria and nucleus, with translation initiation from the upstream start codon allowing for transport to the mitochondria and translation initiation from a downstream start codon allowing for transport to the nucleus. Additionally, alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Recommended Dilutions

WB 1:500 - 1:2000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

Immunogen Information

Gene ID

3980

Swiss Prot

P49916

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 730-1009 of human DNA Ligase III/LIG3 (NP_039269.2).

Synonyms

LIG2; MTDPS20; LIG3alpha; DNA Ligase III/LIG3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using DNA Ligase III/LIG3 Rabbit pAb (A1887) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.

Immunoprecipitation analysis of 200 µg extracts of Jurkat cells using 1 µg DNA Ligase III/LIG3 antibody (A1887). Western blot was performed from the immunoprecipitate using DNA Ligase III/LIG3 antibody (A1887) at a dilution of 1:1000.

