

A19215

Leader in Biomolecular Solutions for Life Science



SARS-CoV-2 neutralizing antibody

Catalog No.: A19215

7 Publications

Basic Information

Observed MW

Calculated MW

Category

Mouse Monoclonal Antibody

Applications

Cross-Reactivity

Background

Recommended Dilutions

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Synonyms

Contact



www.abclonal.com

Product Information


Source

Isotype


Purification

Storage

Validation Data



Immobilized SARS-CoV-2 S1 protein, His Tag (Cat. No. RP01262) at 1 µg/mL (100µL/well) can bind Anti-SARS-CoV-2 Neutralizing Antibody, Human IgG1 (Cat. No. A19215) with a linear range of 0.078-5 ng/mL



Serial dilutions of Anti-SARS-CoV-2 Neutralizing Antibody, Human IgG1 (Cat.No. A19215) was detected by SARS-CoV-2 Inhibitor screening Kit (Cat.No. RK04149) with a half maximal inhibitory concentration (IC50) of 0.496ug/mL.