

A19228

Leader in Biomolecular Solutions for Life Science



LCAT Rabbit mAb

Catalog No.: A19228

Recombinant

Basic Information

Observed MW

62kDa

Calculated MW

50kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

CloneNo number

ARC2383

Background

This gene encodes the extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. Mutations in this gene have been found to cause fish-eye disease as well as LCAT deficiency.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

3931

Swiss Prot

P04180

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LCAT (P04180).

Synonyms

LCAT; lecithin-cholesterol acyltransferase

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of mouse plasma, using LCAT antibody (A19228) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 3s.

Western blot analysis of extracts of rat plasma, using LCAT antibody (A19228) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.

