

A19583

Leader in Biomolecular Solutions for Life Science



Glucocorticoid Receptor Rabbit mAb

Catalog No.: A19583

Recombinant

2 Publications

Basic Information

Observed MW

95kDa

Calculated MW

86kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB,IP,ELISA

Cross-Reactivity

Human,Mouse,Rat

CloneNo number

ARC0062

Background

This gene encodes glucocorticoid receptor, which can function both as a transcription factor that binds to glucocorticoid response elements in the promoters of glucocorticoid responsive genes to activate their transcription, and as a regulator of other transcription factors. This receptor is typically found in the cytoplasm, but upon ligand binding, is transported into the nucleus. It is involved in inflammatory responses, cellular proliferation, and differentiation in target tissues. Mutations in this gene are associated with generalized glucocorticoid resistance. Alternative splicing of this gene results in transcript variants encoding either the same or different isoforms. Additional isoforms resulting from the use of alternate in-frame translation initiation sites have also been described, and shown to be functional, displaying diverse cytoplasm-to-nucleus trafficking patterns and distinct transcriptional activities (PMID:15866175).

Recommended Dilutions

WB 1:500 - 1:2000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

Immunogen Information

Gene ID

2908

Swiss Prot

P04150

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 2-180 of human GR (P04150).

Synonyms

GR; GCR; GRL; GCCR; GCRST; Glucocorticoid Receptor

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using Glucocorticoid Receptor Rabbit mAb (A19583) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.

Immunoprecipitation analysis of 300 µg extracts of U-251MG cells using 3 µg Glucocorticoid Receptor antibody (A19583). Western blot was performed from the immunoprecipitate using Glucocorticoid Receptor antibody (A19583) at a dilution of 1:1000.

