

A19958

Leader in Biomolecular Solutions for Life Science



[KO Validated] ANKRD52 Rabbit pAb

Catalog No.: A19958 **KO** Validated

Basic Information

Observed MW

115kDa

Calculated MW

115kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Putative regulatory subunit of protein phosphatase 6 (PP6) that may be involved in the recognition of phosphoprotein substrates.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

283373

Swiss Prot

Q8NB46

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ANKRD52 (NP_775866.2).

Synonyms

ANKRD33; 52

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from wild type (WT) and ANKRD52 knockout (KO) 293T cells, using [KO Validated] ANKRD52 Rabbit pAb (A19958) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.