

A20016

Leader in Biomolecular Solutions for Life Science



SARS-CoV-2 N Protein Rabbit pAb

Catalog No.: A20016

2 Publications

Basic Information

Observed MW

Calculated MW

46kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,IP,FC,ELISA

Cross-Reactivity

SARS-CoV-2

Background

Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays a fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication. May modulate transforming growth factor-beta signaling by binding host SMAD3.

Recommended Dilutions

WB	1:1000 - 1:5000
ELISA	1:1000 - 1:5000
FC	1:50 - 1:200
IF/ICC	1:50 - 1:200
IP	1:50 - 1:200

Immunogen Information

Gene ID

1489678

Swiss Prot

P59595

Immunogen

Recombinant fusion protein of SARS-CoV-2 Nucleoprotein.

Synonyms

Nucleoprotein; NP; SARS-CoV-2 N Protein

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.