

A20187

Leader in Biomolecular Solutions for Life Science



COL7A1 Rabbit pAb

Catalog No.: A20187

Basic Information

Observed MW

Refer to figures

Calculated MW

295kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes the alpha chain of type VII collagen. The type VII collagen fibril, composed of three identical alpha collagen chains, is restricted to the basement zone beneath stratified squamous epithelia. It functions as an anchoring fibril between the external epithelia and the underlying stroma. Mutations in this gene are associated with all forms of dystrophic epidermolysis bullosa. In the absence of mutations, however, an acquired form of this disease can result from an autoimmune response made to type VII collagen.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

1294

Swiss Prot

Q02388

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 2780-2944 of human COL7A1 (NP_000085.1).

Synonyms

EBD1; EBR1; EBDCT; NDNC8; COL7A1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.