A20202

Leader in Biomolecular Solutions for Life Science

SARS-CoV-2 NSP8 Rabbit pAb

Catalog No.: A20202 2 Publications



Basic Information

Observed MW 25kDa

Calculated MW 141kDa

Category Polyclonal Antibody

Applications WB,ELISA

Cross-Reactivity SARS-CoV-2

Background

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positivesense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ ORF1ab, the largest gene, contains overlapping open reading frames that encode polyproteins PP1ab and PP1a. The polyproteins are cleaved to yield 16 nonstructural proteins, NSP1-16. Production of the longer (PP1ab) or shorter protein (PP1a) depends on a -1 ribosomal frameshifting event. The proteins, based on similarity to other coronaviruses, include the papain-like proteinase protein (NSP3), 3C-like proteinase (NSP5), RNA-dependent RNA polymerase (NSP12, RdRp), helicase (NSP13, HEL), endoRNAse (NSP15), 2'-O-Ribose-Methyltransferase (NSP16) and other nonstructural proteins. SARS-CoV-2 nonstructural proteins are responsible for viral transcription, replication, proteolytic processing, suppression of host immune responses and suppression of host gene expression. The RNAdependent RNA polymerase is a target of antiviral therapies.

Recommended Dilutions

Immunogen Information

WB

1:2000 - 1:10000

Gene ID 43740578 Swiss Prot P0DTC2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of coronavirus NSP8 (YP 009742615.1).

Synonyms

Contact

Ð

al.com Source

www.abclonal.com

Product Information

Purification Affinity purification

Storage


Rabbit

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Isotype

lgG

Validation Data



Western blot analysis of extracts of normal 293T cells and 293T transfected with NSP8 Protein, using SARS-CoV-2 NSP8 Rabbit pAb (A20202) at 1:10000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.