

A20239

Leader in Biomolecular Solutions for Life Science



PLXND1 Rabbit pAb

Catalog No.: A20239

Basic Information

Observed MW

Refer to figures

Calculated MW

212kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables protein domain specific binding activity. Predicted to be involved in several processes, including endothelial cell migration; nervous system development; and regulation of angiogenesis. Predicted to act upstream of or within several processes, including circulatory system development; dichotomous subdivision of terminal units involved in salivary gland branching; and positive regulation of protein binding activity. Located in lamellipodium.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23129

Swiss Prot

Q9Y4D7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human PLXND1 (NP_055918.2).

Synonyms

CHTD9; PLEXD1; PLXND1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.