

A20243

Leader in Biomolecular Solutions for Life Science



RRAS Rabbit pAb

Catalog No.: A20243

Basic Information

Observed MW

Refer to figures

Calculated MW

23kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is a small GTPase involved in diverse processes including angiogenesis, vascular homeostasis and regeneration, cell adhesion, and neuronal axon guidance. Mutations in this gene are found in many invasive cancers.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

6237

Swiss Prot

P10301

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human RRAS (NP_006261.1).

Synonyms

R-Ras; RRAS

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.