

A20297

Leader in Biomolecular Solutions for Life Science



TYRP1 Rabbit pAb

Catalog No.: A20297

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

61kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a melanosomal enzyme that belongs to the tyrosinase family and plays an important role in the melanin biosynthetic pathway. Defects in this gene are the cause of rufous oculocutaneous albinism and oculocutaneous albinism type III.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

7306

Swiss Prot

P17643

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 350-450 of human TYRP1 (NP_000541.1).

Synonyms

TRP; CAS2; CATB; GP75; OCA3; TRP1; TYRP; b-PROTEIN; TYRP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.