

A20300

Leader in Biomolecular Solutions for Life Science



ATG101 Rabbit pAb

Catalog No.: A20300

Basic Information

Observed MW

Refer to figures

Calculated MW

25kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Involved in autophagosome assembly. Located in phagophore assembly site.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

60673

Swiss Prot

Q9BSB4

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 100-200 of human ATG101 (NP_001092143.1).

Synonyms

C12orf44; ATG101

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.