

A20354

Leader in Biomolecular Solutions for Life Science



TNFAIP8L1 Rabbit pAb

Catalog No.: A20354

Basic Information

Observed MW

Refer to figures

Calculated MW

21kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables identical protein binding activity. Predicted to be involved in negative regulation of TOR signaling. Predicted to be active in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

126282

Swiss Prot

Q8WVP5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-186 of human TNFAIP8L1 (NP_689575.2).

Synonyms

TIPE1; TNFAIP8L1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.