

A20358

Leader in Biomolecular Solutions for Life Science



CCDC85C Rabbit pAb

Catalog No.: A20358

Basic Information

Observed MW

Refer to figures

Calculated MW

45kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in cerebral cortex development. Located in adherens junction.

Recommended Dilutions

Elisa 1:500 - 1:2000

WB 1:500 - 1:2000

Immunogen Information

Gene ID

317762

Swiss Prot

A6NKD9

Immunogen

Synthetic peptide of human CCDC85C

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% sodium azide, 40% glycerol, pH 7.4.