

A20359

Leader in Biomolecular Solutions for Life Science



SHROOM1 Rabbit pAb

Catalog No.: A20359

Basic Information

Observed MW

100-130kDa

Calculated MW

91kDa

Category

Polyclonal Antibody

Applications

WB,IP,ELISA

Cross-Reactivity

Human

Background

SHROOM family members play diverse roles in the development of the nervous system and other tissues (Hagens et al., 2006 [PubMed 16615870]).

Recommended Dilutions

Elisa	1:500 - 1:2000
WB	1:500 - 1:3000
IP	1:500 - 1:1000

Immunogen Information

Gene ID

134549

Swiss Prot

Immunogen

Recombinant protein of human SHROOM1.

Synonyms

APXL2; SHROOM1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.