

A20361

Leader in Biomolecular Solutions for Life Science



## RASSF7 Rabbit pAb

Catalog No.: A20361

### Basic Information

**Observed MW**

40kDa

**Calculated MW**

40kDa

**Category**

Mouse Monoclonal Antibody

**Applications**

WB,ELISA

**Cross-Reactivity**

Human,Rat

### Background

Predicted to be involved in apoptotic process; regulation of microtubule cytoskeleton organization; and signal transduction. Located in centriolar satellite.

### Recommended Dilutions

**Elisa** 1:500 - 1:2000

**WB** 1:500 - 1:1000

### Immunogen Information

**Gene ID**

8045

**Swiss Prot**

**Immunogen**

Recombinant protein of human RASSF7.

**Synonyms**

HRC1; FAP88; HRAS1; CFAP88; C11orf13; RASSF7

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.