

A20362

Leader in Biomolecular Solutions for Life Science



SINHCAF Rabbit pAb

Catalog No.: A20362

Basic Information

Observed MW

Refer to figures

Calculated MW

25kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Involved in negative regulation of cell migration. Part of Sin3 complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

58516

Swiss Prot

Q9NP50

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SINHCAF (NP_001129283.1).

Synonyms

L4; TERA; FAM60A; C12orf14; SINHCAF

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.