

# ZNF318 Rabbit pAb

Catalog No.: A20374

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

251kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

Predicted to enable protein heterodimerization activity and protein homodimerization activity. Predicted to be involved in negative regulation of transcription, DNA-templated and positive regulation of transcription, DNA-templated. Located in cytosol and nucleoplasm.

## Recommended Dilutions

**WB** 1:200 - 1:2000

## Immunogen Information

### Gene ID

24149

### Swiss Prot

Q5VUA4

### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ZNF318 (NP\_055160.2).

### Synonyms

TZF; ZFP318; HRIHFB2436; ZNF318

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.