

A20376

Leader in Biomolecular Solutions for Life Science



## EXOSC9 Rabbit pAb

Catalog No.: A20376

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

49kDa

#### Category

Polyclonal Antibody

#### Applications

WB, IF/ICC

#### Cross-Reactivity

Human

### Background

This gene encodes a component of the human exosome, a exoribonuclease complex which processes and degrades RNA in the nucleus and cytoplasm. This component may play a role in mRNA degradation and the polymyositis/scleroderma autoantigen complex. Alternative splicing results in multiple transcript variants.

### Recommended Dilutions

WB	1:1000 - 1:6000
IF/ICC	1:10 - 1:100

### Immunogen Information

#### Gene ID

5393

#### Swiss Prot

Q06265

#### Immunogen

Recombinant protein of human EXOSC9.

#### Synonyms

p5; p6; PCH1D; RRP45; PMSCL1; Rrp45p; PM/Scl-75; EXOSC9

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% sodium azide, 50% glycerol, pH 7.3.