

A20384

Leader in Biomolecular Solutions for Life Science



Fam221b Rabbit pAb

Catalog No.: A20384

Basic Information

Observed MW

55kDa

Calculated MW

54kDa

Category

Polyclonal Antibody

Applications

WB,IHC-P,ELISA

Cross-Reactivity

Mouse,Rat

Background

Orthologous to human FAM221B (family with sequence similarity 221 member B).

Recommended Dilutions

WB 1:500 - 1:1000

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

242408

Swiss Prot

Q8C627

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 280-487 of mouse Fam221b (NP_780726.1).

Synonyms

4930412F15Rik; Fam221b

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using Fam221b Rabbit pAb (A20384) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.



Immunohistochemistry analysis of Fam221b in paraffin-embedded rat testis using Fam221b Rabbit pAb (A20384) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.