

IFIT5 Rabbit pAb

Catalog No.: A20391

Basic Information

Observed MW

Refer to figures

Calculated MW

56kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables nucleic acid binding activity. Involved in defense response to virus; negative regulation of viral genome replication; and positive regulation of I-kappaB kinase/NF-kappaB signaling. Located in actin cytoskeleton; apical part of cell; and ruffle membrane. Colocalizes with IkappaB kinase complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

24138

Swiss Prot

Q13325

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 183-482 of human IFIT5 (NP_036552.1).

Synonyms

P58; RI58; ISG58; IFIT5

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.