

A20408

Leader in Biomolecular Solutions for Life Science



MRM1 Rabbit pAb

Catalog No.: A20408

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables rRNA (guanosine-2'-O-)-methyltransferase activity. Predicted to be involved in rRNA 2'-O-methylation. Located in mitochondrion.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

79922

Swiss Prot

Q6IN84

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human MRM1 (NP_079140.2).

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.