

A20411

Leader in Biomolecular Solutions for Life Science



RBM23 Rabbit pAb

Catalog No.: A20411

Basic Information

Observed MW

Refer to figures

Calculated MW

49kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the U2AF-like family of RNA binding proteins. This protein interacts with some steroid nuclear receptors, localizes to the promoter of a steroid-responsive gene, and increases transcription of steroid-responsive transcriptional reporters in a hormone-dependent manner. It is also implicated in the steroid receptor-dependent regulation of alternative splicing. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

55147

Swiss Prot

Q86U06

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human RBM23 (NP_060577.3).

Synonyms

PP239; RNPC4; CAPERbeta; RBM23

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.