

A20417

Leader in Biomolecular Solutions for Life Science



RGL2 Rabbit pAb

Catalog No.: A20417

Basic Information

Observed MW

Refer to figures

Calculated MW

84kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable guanyl-nucleotide exchange factor activity. Predicted to be involved in regulation of catalytic activity and small GTPase mediated signal transduction. Predicted to act upstream of or within negative regulation of cardiac muscle cell apoptotic process; positive regulation of phosphatidylinositol 3-kinase signaling; and regulation of Ral protein signal transduction. Predicted to be located in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5863

Swiss Prot

O15211

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RGL2 (NP_004752.1).

Synonyms

KE1.5; RAB2L; HKE1.5; RGL2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.