

A20422

Leader in Biomolecular Solutions for Life Science



THAP11 Rabbit pAb

Catalog No.: A20422

Basic Information

Observed MW

Refer to figures

Calculated MW

34kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene contains a THAP domain, which is a conserved DNA-binding domain that has striking similarity to the site-specific DNA-binding domain (DBD) of Drosophila P element transposases.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

57215

Swiss Prot

Q96EK4

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 200-300 of human THAP11 (NP_065190.2).

Synonyms

RONIN; CTG-B43a; CTG-B45d; HRIHFB2206; THAP11

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.