

A20429

Leader in Biomolecular Solutions for Life Science



# GATAD1 Rabbit pAb

Catalog No.: A20429

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

29kDa

### Category

Polyclonal Antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human

## Background

The protein encoded by this gene contains a zinc finger at the N-terminus, and is thought to bind to a histone modification site that regulates gene expression. Mutations in this gene have been associated with autosomal recessive dilated cardiomyopathy. Alternatively spliced transcript variants have been found for this gene.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

57798

### Swiss Prot

Q8WUU5

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-269 of human GATAD1 (NP\_066990.3).

### Synonyms

ODAG; CMD2B; RG083M05.2; GATAD1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.