

PRSS35 Rabbit pAb

Catalog No.: A20430

Basic Information

Observed MW

Refer to figures

Calculated MW

47kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to be located in extracellular region.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

167681

Swiss Prot

Q8N3Z0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 17-270 of human PRSS35 (NP_699193.2).

Synonyms

C6orf158; dj223E3.1; PRSS35

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.