

A20434

Leader in Biomolecular Solutions for Life Science



AMER1 Rabbit pAb

Catalog No.: A20434

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

124kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene upregulates transcriptional activation by the Wilms tumor protein and interacts with many other proteins, including CTNNB1, APC, AXIN1, and AXIN2. Defects in this gene are a cause of osteopathia striata with cranial sclerosis (OSCS).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

139285

Swiss Prot

Q5JTC6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 550-660 of human AMER1 (NP_689637.3).

Synonyms

WTX; OSCS; FAM123B; AMER1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.